
| PRODUCT NAME | 4 MOTOR DIALYSIS SEAT |
|-----------------------------|--|
| ORIGINAL PRODUCT CODE (P/N) | NHS 900 TERRA - 4 MOTORS DIALYSIS CHAIR (LINAK COLUMN) |

PRODUCT PICTURE



TECHNICAL DRAWING



| | ١٧S | CF | AT |
|--|-----|----|----|
| | | | |

| DIALYSIS SEAT | | | | |
|--|--|--|--|--|
| Type(Manuel-Motor) | DIALYSIS SEAT WITH MOTOR | | | |
| Height Adjustment | MIN : 590 , MAX : 990 | | | |
| Back Leght | Lenght : 790 EN : 550 | | | |
| Seating Lenght | Lenght : 480 EN : 600 | | | |
| Foot Long | Lenght: 570. Consists of 2 Parts eith foot stanting side. | | | |
| Back Angle Adjustment | It has.By help of Motor. 0-70° | | | |
| Trendelenburg-reverse Trendelenburg angle adjustment | It has .by help of motor 0- 13° | | | |
| Moveable Sides mechanism (Remote Control System) | It has remote control. Seat height adjustment, Back movement, Foot movement, TRENDELENBURG and Foot extention s made by Remote Control | | | |
| Foot Angle | Help of Motor from 0- to 40° | | | |
| Accessories | IV POL MADE BY STAINESS STEEL, SPECIAL VINLEKS MATERIAL COVERED PILLOW. ON THE SEATING SITE 5 CM VISCO PILLOWWAS UDES IN ORDER TO MAKE YOU FEEL COMFORTABLE. AYRICA EST SETTLE DOWM MECHANISM FOOT PUTTING SIDE CAN EXTEND 20 CM ACCORDING TO PATIENT HEIGHT.ARMS CAN BE ADJUSTED TO RIGHT, LEFT AND UP DOWN. ARMS CAN BE ADJUSTED FOR PATIENT COMFORT.MOVEABLE PARTS ARE MADE BY STAINLESS STEEL. | | | |
| Production Accessories | UP & DOWN ADJUSTMENT | | | |
| Armlet Width | Width: 200 Height: 680 VINLEKS COVERED MATERIAL, ABS METAL SHEET AND ALUMINUM CONNECTION | | | |
| lv pol | Single hook metal thickness is 1,5 cm | | | |
| Hand Handle | Back side has 25mm thickness metal handle | | | |
| Carrying Chasis Material | 30 x 40 x 2 mm ,30 x 50 x 2,20 x 20 x 2 and 20 x 30 x 2 mm BOX PROFILE , Ø 25 x 2 PIPE PROFILE ve 5 mm, 6 mm and 8 mm METAL SHEETCHASIS IS MADE BY CARBON STEEL. ELECTROSTATIC POWDER PAINT IS USED AGAINST RUST. | | | |
| Back Setiang And Resting Side Material | Seating side covered antibacterial material. | | | |
| SAFE LOADING | 225 KG | | | |

| Castors | DIAMETER 125 mm / IMPORTED CASTORS / 4 UNIT | | |
|----------------------------------|---|--|--|
| Motor Number and Efficiency Help | 4 MOTOR. Height Adjustment (2000 N), Back Adjustment (3500 N), Foot Adjustment (6000 N), TRENDELENBURG (6000 N) | | |
| Other Technical Detail | Patient Sliping prevention by help of motor ,4 unit central locking pedal, it has acumulator.Compact armlets according to desire. Digital Scale can be Put according to desire in order to observe paint weight | | |